

EARTH 48 SINGLE

Aluminum Body

LED



Number of LED
48

Intensity
160,000 Lux

Color Temp
3800 K - 5200 K

Power
48 Watt

Diameter Of Light
500 mm

Technical Data

Number of LED	48
Intensity	160,000
Spot Diameter	100 - 150 mm
Diameter of Light	500 mm
Color Temperature	3800 - 5200K
Depth illumination	10-15 cm
Color Rendering Index	95 Ra
Focusing	Adjustable
LED Average Life	> 50000 Hrs.
Intensity Control	10 -100 Digital
Power Supply	220V/50Hz AC
Power	48Watt
Special Function	5 memory Function with Remote
Optional	Battery Backup 360 degree structure



Multi functional Controller
With Remote

EARTH 48 MOBILE

Aluminum Body

LED



Number of LED
48

Intensity
160,000 Lux

Color Temp
3800 K - 5200 K

Power
48 Watt

Diameter Of Light
500 mm



Multi functional Controller
With Remote



Optional Battery Backup

Technical Data

Number of LED	48	Focusing	Adjustable
Intensity	160, 000 Lux	LED Average Life	> 50000 Hrs.
Spot Diameter	100 - 150 mm	Intensity Control	10 -100 Digital
Color Temperature	3800 - 5200K	Power Supply	220V/50Hz AC
Depth illumination	10-15 cm	Power	48 Watt
Color Rendering Index	95 Ra	Special Function	5 memory Function with Remote
Diameter of Light	500 mm		

