

GLANZA 500 SINGLE

LED

Number of LED
48



Intensity
160,000 Lux



Color Temp
3800 K - 5200 K



Power
48 Watt



Diameter Of Light
500 mm

Technical Data

Number of LED	48
Intensity	160,000
Spot Diameter	100 - 150 mm
Diameter of Light	500 mm
Color Temperature	3800 - 5200K
Depth illumination	10-15 cm
Color Rendering Index	95 Ra
Focusing	Adjustable
LED Average Life	> 50000 Hrs.
Intensity Control	10 -100 Digital
Power Supply	220V/50Hz AC
Power	48Watt
Special Function	5 memory Function
	with Remote
Optional	Battery Backup
	360 degree structure
Body Type	Aluminium



Multi functional Controller
With Remote

GLANZA 500 MOBILE

LED

Number of LED
48



Intensity
160,000 Lux



Color Temp
3800 K - 5200 K



Power
48 Watt



Diameter Of Light
500 mm

Technical Data

Number of LED	48	Focusing	Adjustable
Intensity	160, 000 Lux	LED Average Life	> 50000 Hrs.
Spot Diameter	100 - 150 mm	Intensity Control	10 -100 Digital
Color Temperature	3800 - 5200K	Power Supply	220V/50Hz AC
Depth illumination	10-15 cm	Power	48 Watt
Color Rendering Index	95 Ra	Special Function	5 memory Function
Diameter of Light	500 mm		with Remote
		Body Type	Aluminium



Multi functional Controller
With Remote



Optional Battery Backup

